

REOPRO BOLUS & INFUSION DOSING CHART

(For patients with acute coronary syndromes presenting immediately to CCL)

Based on Bolus dose of 0.25 mg/kg and

Infusion dose of 0.125 mcg/kg/min (max of 10 mcg/min)

All infusions are run @ 21ml/hr for 12 hours

Wt in kg	Bolus in ml	ml to add to 250ml		Wt in kg	Bolus in ml	ml to add to 250ml		Wt in kg	Bolus in ml	ml to add to 250ml
45	5.6	2.03		80	10	3.6		115	14.4	3.6
46	5.8	2.07		81	10.1	3.6		116	14.5	3.6
47	5.9	2.12		82	10.2	3.6		117	14.6	3.6
48	6.0	2.16		83	10.4	3.6		118	14.8	3.6
49	6.1	2.21		84	10.5	3.6		119	14.9	3.6
50	6.2	2.25		85	10.6	3.6		120	15.0	3.6
51	6.4	2.3		86	10.8	3.6		121	15.1	3.6
52	6.5	2.34		87	10.9	3.6		122	15.3	3.6
53	6.6	2.39		88	11.0	3.6		123	15.4	3.6
54	6.8	2.43		89	11.1	3.6		124	15.5	3.6
55	6.9	2.48		90	11.2	3.6		125	15.6	3.6
56	7.0	2.52		91	11.4	3.6		126	15.8	3.6
57	7.1	2.57		92	11.5	3.6		127	15.9	3.6
58	7.2	2.61		93	11.6	3.6		128	16.0	3.6
59	7.4	2.66		94	11.8	3.6		129	16.1	3.6
60	7.5	2.7		95	11.9	3.6		130	16.3	3.6
61	7.6	2.75		96	12.0	3.6		131	16.4	3.6
62	7.8	2.79		97	12.1	3.6		132	16.5	3.6
63	7.9	2.84		98	12.3	3.6		133	16.6	3.6
64	8.0	2.88		99	12.4	3.6		134	16.8	3.6
65	8.1	2.93		100	12.5	3.6		135	16.9	3.6
66	8.2	2.97		101	12.6	3.6		136	17.0	3.6
67	8.4	3.02		102	12.8	3.6		137	17.1	3.6
68	8.5	3.06		103	12.9	3.6		138	17.3	3.6
69	8.6	3.11		104	13.0	3.6		139	17.4	3.6
70	8.8	3.15		105	13.1	3.6		140	17.5	3.6
71	8.9	3.2		106	13.3	3.6		141	17.6	3.6
72	9.0	3.24		107	13.4	3.6		142	17.8	3.6
73	9.1	3.29		108	13.5	3.6		143	17.9	3.6
74	9.2	3.33		109	13.6	3.6		144	18.0	3.6
75	9.4	3.38		110	13.8	3.6		145	18.1	3.6
76	9.5	3.42		111	13.9	3.6		146	18.3	3.6
77	9.6	3.47		112	14	3.6		147	18.4	3.6
78	9.8	3.51		113	14.1	3.6		148	18.5	3.6
79	9.9	3.56		114	14.3	3.6				

- **Withdraw amount for bolus & infusion from Reopro vial & filter with 0.22 µm low protein-binding filter. Mix infusion in 250 ml NS. Rate = 21 ml/hr x 12 hours only.**
- **DO NOT SHAKE!!!!!!!**